



ASSURE-100

Rapid COVID-19 Test



FDA EUA Authorized



ASSURE-100 Rapid COVID-19 Test - SIMPLE, FAST, ACCURATE

The ASSURE-100 Rapid COVID-19 Test conveniently delivers an accurate and reliable COVID-19 diagnosis in 20 minutes. It is a point-of-care tool that detects SARS-CoV-2 antigens directly from shallow nasal swab specimens and is visually read, not requiring any instrumentation.

SIMPLE – Swab, Mix, Pour, interpret Results.

FAST – Accurate and reliable results within 20 minutes, allowing for COVID-19 testing of patients within the first eight days of symptoms.

ACCURATE – ASSURE-100 has a sensitivity (PPA) of 89% and specificity (NPA) of 100% and can detect all major variants of COVID-19.



About Oceanit: Founded in 1985 in Honolulu, Hawaii, Oceanit has earned a world-class reputation for breakthrough science & technology, sustainable engineering, and disruptive innovation. Oceanit Foundry LLC is a subsidiary company of Oceanit Labs, Inc.

This product has not been FDA cleared or approved, but has been authorized by FDA under an EUA for use by authorized laboratories; this product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.